

CYP2E1 Antibody, Biotin conjugated



PACO26104

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Metabolizes several precarcinogens, drugs, and solvents to reactive metabolites. Inactivates a number of drugs and xenobiotics and also bioactivates many xenobiotic substrates to their hepatotoxic or carcinogenic forms.

Gene ID:

CYP2E1

Uniprot

P05181

Synonyms:

Cytochrome P450 2E1 (EC 1.14.13. -) (4-nitrophenol 2-hydroxylase) (EC 1.14.13. n7) (CYP11E1) (Cytochrome P450-J), CYP2E1, CYP2E

Immunogen:

Recombinant Human Cytochrome P450 2E1 protein (49-288AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A