

CYP3A7 Antibody, HRP conjugated



PACO26106

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway. It oxidizes a variety of structurally unrelated compounds, including steroids, fatty acid, , and xenobiotics.

Gene ID:

CYP3A7

Uniprot

P24462

Synonyms:

Cytochrome P450 3A7 (EC 1.14.14.1) (CYP11A7) (Cytochrome P450-HFLA), CYP3A7

Immunogen:

Recombinant Human Cytochrome P450 3A7 protein (1-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A