

RPS25 Antibody



PACO26109

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

RPS25

Uniprot

P62851

Synonyms:

40S ribosomal protein S25 (Small ribosomal subunit protein eS25), RPS25

Immunogen:

Recombinant Human 40S ribosomal protein S25 protein (1-125AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A