

RPS25 Antibody, Biotin conjugated



PACO26112

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	RPS25
Source:	Uniprot
Rabbit	P62851
Isotype:	Synonyms:
IgG	40S ribosomal protein S25 (Small ribosomal subunit protein eS25), RPS25
Applications:	Immunogen:
ELISA	Recombinant Human 40S ribosomal protein S25 protein (1-125AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A