

DAOA Antibody, FITC conjugated



PACO26127

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	DAOA
Source:	Uniprot
Rabbit	P59103
Isotype:	Synonyms:
IgG	D-amino acid, oxidase activator (Protein G72), DAOA, G72
Applications:	Immunogen:
ELISA	Recombinant Human D-amino acid, oxidase activator protein (1-153AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A