

Dcn Antibody, FITC conjugated



PACO26158

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May affect the rate of fibrils formation.

Gene ID:

Dcn

Uniprot

P28654

Synonyms:

Decorin (Bone proteoglycan II) (PG-S2) (PG40), Dcn

Immunogen:

Recombinant Mouse Decorin protein (35-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A