

# DCTN3 Antibody, FITC conjugated



PACO26167

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Together with dynein may be involved in spindle assembly and cytokinesis.

**Gene ID:**

DCTN3

**Uniprot**

O75935

**Synonyms:**

Dynactin subunit 3 (Dynactin complex subunit 22 kDa subunit) (p22), DCTN3, DCTN22

**Immunogen:**

Recombinant Human Dynactin subunit 3 protein (2-176AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A