

# DGKG Antibody, HRP conjugated



PACO26226

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Reverses the normal flow of glycerolipid biosynthesis by phosphorylating diacylglycerol back to phosphatidic acid, .

**Gene ID:**

DGKG

**Uniprot**

P49619

**Synonyms:**

Diacylglycerol kinase gamma (DAG kinase gamma) (EC 2.7.1.107) (Diglyceride kinase gamma) (DGK-gamma), DGKG, DAGK3

**Immunogen:**

Recombinant Human Diacylglycerol kinase &gamma; protein (1-255AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A