

DHRS3 Antibody, HRP conjugated



PACO26238

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the reduction of all-trans-retinal to all-trans-retinol in the presence of NADPH.

Gene ID:

DHRS3

Uniprot

O75911

Synonyms:

Short-chain dehydrogenase/reductase 3 (EC 1.1.1.300) (DD83.1) (Retinal short-chain dehydrogenase/reductase 1) (retSDR1) (Retinol dehydrogenase 17) (Short chain dehydrogenase/reductase family 16C member 1), DHRS3, RDH17 SDR16C1

Immunogen:

Recombinant Human Short-chain dehydrogenase/reductase 3 protein (30-169AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A