

DNAJB2 Antibody, HRP conjugated



PACO26266

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

DNAJB2

Uniprot

P25686

Synonyms:

DnaJ homolog subfamily B member 2 (Heat shock 40 kDa protein 3) (Heat shock protein J1) (HSJ-1), DNAJB2, HSJ1 HSPF3

Immunogen:

Recombinant Human DnaJ homolog subfamily B member 2 protein (2-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A