

DNAL4 Antibody, HRP conjugated



PACO26278

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Force generating protein of respiratory cilia. Produces force towards the minus ends of microtubules. Dynein has ATPase activity.

Gene ID:

DNAL4

Uniprot

O96015

Synonyms:

Dynein light chain 4, axonemal, DNAL4

Immunogen:

Recombinant Human Dynein light chain 4, axonemal protein (1-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A