

# DNALI1 Antibody, HRP conjugated



PACO26282

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a dynamic role in flagellar motility.

**Gene ID:**

DNALI1

**Uniprot**

O14645

**Synonyms:**

Axonemal dynein light intermediate polypeptide 1 (Inner dynein arm light chain, axonemal) (hp28), DNALI1

**Immunogen:**

Recombinant Human Axonemal dynein light intermediate polypeptide 1 protein (1-258AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A