

B4GALT3 Antibody, FITC conjugated



PACO26287

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

B4GALT3

Uniprot

O60512

Synonyms:

Beta-1,4-galactosyltransferase 3 (Beta-1,4-GalTase 3) (Beta4Gal-T3) (b4Gal-T3) (EC 2.4.1.-) (UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 3) (UDP-galactose: beta-N-acetylglucosamine beta-1,4-galactosyltransferase 3) [Includes: N-acetylglucosamine synthase (EC 2.4.1.90) (Nal synthase); Beta-N-acetylglucosaminylglycopeptide beta-1,4-galactosyltransferase (EC 2.4.1.38); Beta-N-acetylglucosaminyl-glycolipid beta-1,4-galactosyltransferase (EC 2.4.1.-)], B4GALT3

Immunogen:

Recombinant Human Beta-1,4-galactosyltransferase 3 protein (34-393AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A