

DNTTIP1 Antibody, FITC conjugated



PACO26307

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

DNTTIP1

Uniprot

Q9H147

Synonyms:

Deoxynucleotidyltransferase terminal-interacting protein 1 (Terminal deoxynucleotidyltransferase-interacting factor 1) (TdIF1) (TdT-interacting factor 1), DNTTIP1, C20orf167 TDIF1

Immunogen:

Recombinant Human Deoxynucleotidyltransferase terminal-interacting protein 1 protein (1-329AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A