

# DNTTIP1 Antibody, Biotin conjugated



PACO26308

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

DNTTIP1

**Uniprot**

Q9H147

**Synonyms:**

Deoxynucleotidyltransferase terminal-interacting protein 1 (Terminal deoxynucleotidyltransferase-interacting factor 1) (TdIF1) (TdT-interacting factor 1), DNTTIP1, C20orf167 TDIF1

**Immunogen:**

Recombinant Human Deoxynucleotidyltransferase terminal-interacting protein 1 protein (1-329AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A