

# DOK5 Antibody, HRP conjugated



PACO26310

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

DOK proteins are enzymatically inert adaptor or scaffolding proteins. They provide a docking platform for the assembly of multimolecular signaling complexes. DOK5 functions in RET-mediated neurite outgrowth and plays a positive role in activation of the MAP kinase pathway. Putative link with downstream effectors of RET in neuronal differentiation.

**Gene ID:**

DOK5

**Uniprot**

Q9P104

**Synonyms:**

Docking protein 5 (Downstream of tyrosine kinase 5) (Insulin receptor substrate 6) (IRS-6) (IRS6), DOK5, C20orf180

**Immunogen:**

Recombinant Human Docking protein 5 protein (1-198AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A