

# DRD5 Antibody, Biotin conjugated



PACO26324

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Dopamine receptor whose activity is mediated by G proteins which activate adenylyl cyclase.

**Gene ID:**

DRD5

**Uniprot**

P21918

**Synonyms:**

D(1B) dopamine receptor (D(5) dopamine receptor) (D1beta dopamine receptor) (Dopamine D5 receptor), DRD5, DRD1B DRD1L2

**Immunogen:**

Recombinant Human D(1B) dopamine receptor protein (361-477AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A