

DSCR4 Antibody



PACO26329

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

DSCR4

Uniprot

P56555

Synonyms:

Down syndrome critical region protein 4 (Down syndrome critical region protein B), DSCR4, DCRB DSCRB

Immunogen:

Recombinant Human Down syndrome critical region protein 4 protein (1-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A