

# DSP Antibody, FITC conjugated



PACO26347

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Major high molecular weight protein of desmosomes. Involved in the organization of the desmosomal cadherin-plakoglobin complexes into discrete plasma membrane domains and in the anchoring of intermediate filaments to the desmosomes.

**Gene ID:**

DSP

**Uniprot**

P15924

**Synonyms:**

Desmoplakin (DP) (250/210 kDa paraneoplastic pemphigus antigen), DSP

**Immunogen:**

Recombinant Human Desmoplakin protein (78-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A