

DSP Antibody, Biotin conjugated



PACO26348

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Major high molecular weight protein of desmosomes. Involved in the organization of the desmosomal cadherin-plakoglobin complexes into discrete plasma membrane domains and in the anchoring of intermediate filaments to the desmosomes.

Gene ID:

DSP

Uniprot

P15924

Synonyms:

Desmoplakin (DP) (250/210 kDa paraneoplastic pemphigus antigen), DSP

Immunogen:

Recombinant Human Desmoplakin protein (78-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A