

PACO26357

Product Information

Size:

50ug

Reactivity:

Bacillus anthracis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Phosphorylation of dTMP to form dTDP in both de novo and salvage pathways of dTTP synthesis By similarity.

Gene ID:

tmk

Uniprot

C3LJ02

Synonyms:

Thymidylate kinase (EC 2.7.4.9) (dTMP kinase), tmk

Immunogen:

Recombinant Bacillus anthracis Thymidylate kinase protein (1-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A