

Efna5 Antibody, Biotin conjugated



PACO26404

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Efna5

Uniprot

P97605

Synonyms:

Ephrin-A5 (AL-1) (EPH-related receptor tyrosine kinase ligand 7) (LERK-7), Efna5, Eplg7 Lerk7

Immunogen:

Recombinant Rat Ephrin-A5 protein (21-203AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A