

EGLN2 Antibody, Biotin conjugated



PACO26412

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EGLN2

Uniprot

Q96KS0

Synonyms:

Egl nine homolog 2 (EC 1.14.11.29) (Estrogen-induced tag 6) (HPH-3) (Hypoxia-inducible factor prolyl hydroxylase 1) (HIF-PH1) (HIF-prolyl hydroxylase 1) (HPH-1) (Prolyl hydroxylase domain-containing protein 1) (PHD1), EGLN2, EIT6

Immunogen:

Recombinant Human Egl nine homolog 2 protein (283-407AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A