

# EIF4EBP2 Antibody, HRP conjugated



PACO26422

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Regulates eIF4E activity by preventing its assembly into the eIF4F complex. Mediates the regulation of protein translation by hormones, growth factors and other stimuli that signal through the MAP kinase pathway.

**Gene ID:**

EIF4EBP2

**Uniprot**

Q13542

**Synonyms:**

Eukaryotic translation initiation factor 4E-binding protein 2 (4E-BP2) (eIF4E-binding protein 2), EIF4EBP2

**Immunogen:**

Recombinant Human Eukaryotic translation initiation factor 4E-binding protein 2 protein (1-120AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A