

EIF4EBP2 Antibody, Biotin conjugated



PACO26424

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulates eIF4E activity by preventing its assembly into the eIF4F complex. Mediates the regulation of protein translation by hormones, growth factors and other stimuli that signal through the MAP kinase pathway.

Gene ID:

EIF4EBP2

Uniprot

Q13542

Synonyms:

Eukaryotic translation initiation factor 4E-binding protein 2 (4E-BP2) (eIF4E-binding protein 2), EIF4EBP2

Immunogen:

Recombinant Human Eukaryotic translation initiation factor 4E-binding protein 2 protein (1-120AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A