

ELAVL2 Antibody, Biotin conjugated



PACO26452

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds RNA. Seems to recognize a GAAA motif.

Gene ID:

ELAVL2

Uniprot

Q12926

Synonyms:

ELAV-like protein 2 (ELAV-like neuronal protein 1) (Hu-antigen B) (HuB) (Nervous system-specific RNA-binding protein Hel-N1), ELAVL2, HUB

Immunogen:

Recombinant Human ELAV-like protein 2 protein (4-356AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A