

ELAVL4 Antibody, Biotin conjugated



PACO26456

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in neuron-specific RNA processing. Protects CDKN1A mRNA from decay by binding to its 3'-UTR. Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF and FOS mRNA.

Gene ID:

ELAVL4

Uniprot

P26378

Synonyms:

ELAV-like protein 4 (Hu-antigen D) (HuD) (Paraneoplastic encephalomyelitis antigen HuD), ELAVL4, HUD PNEM

Immunogen:

Recombinant Human ELAV-like protein 4 protein (1-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A