

EMP1 Antibody



PACO26465

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

integral component of membrane, membrane, plasma membrane, bleb assembly, cell death, cell proliferation, epidermis development, multicellular organismal development.

Gene ID:

EMP1

Uniprot

P54849

Synonyms:

Epithelial membrane protein 1 (EMP-1) (CL-20) (Protein B4B) (Tumor-associated membrane protein), EMP1, B4B TMP


Immunogen:

Recombinant Human Epithelial membrane protein 1 protein (22-66AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO26465 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO26465 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).