

ENTPD3 Antibody, HRP conjugated



PACO26490

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has a threefold preference for the hydrolysis of ATP over ADP.

Gene ID:

ENTPD3

Uniprot

O75355

Synonyms:

Ectonucleoside triphosphate diphosphohydrolase 3 (NTPDase 3) (EC 3.6.1.5) (CD39 antigen-like 3) (Ecto-ATP diphosphohydrolase 3) (Ecto-ATPDase 3) (Ecto-ATPase 3) (Ecto-apyrase 3) (HB6), ENTPD3, CD39L3

Immunogen:

Recombinant Human Ectonucleoside triphosphate diphosphohydrolase 3 protein (60-338AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A