

# ENTPD3 Antibody, Biotin conjugated



PACO26492

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Has a threefold preference for the hydrolysis of ATP over ADP.

**Gene ID:**

ENTPD3

**Uniprot**

O75355

**Synonyms:**

Ectonucleoside triphosphate diphosphohydrolase 3 (NTPDase 3) (EC 3.6.1.5) (CD39 antigen-like 3) (Ecto-ATP diphosphohydrolase 3) (Ecto-ATPDase 3) (Ecto-ATPase 3) (Ecto-apyrase 3) (HB6), ENTPD3, CD39L3

**Immunogen:**

Recombinant Human Ectonucleoside triphosphate diphosphohydrolase 3 protein (60-338AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A