

# ENTPD7 Antibody, Biotin conjugated



PACO26496

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Preferentially hydrolyzes nucleoside 5'-triphosphates. The order of activity with respect to possible substrates is UTP > GTP > CTP.

**Gene ID:**

ENTPD7

**Uniprot**

Q9NQZ7

**Synonyms:**

Ectonucleoside triphosphate diphosphohydrolase 7 (NTPDase 7) (EC 3.6.1. -)  
(Lysosomal apyrase-like protein 1), ENTPD7, LALP1

**Immunogen:**

Recombinant Human Ectonucleoside triphosphate diphosphohydrolase 7 protein (51-350AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A