

# EPO Antibody, Biotin conjugated



PACO26520

---

## Product Information

**Size:**

50ug

**Reactivity:**

Macaca fascicularis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Erythropoietin is the principal hormone involved in the regulation of erythrocyte differentiation and the maintenance of a physiological level of circulating erythrocyte mass.

**Gene ID:**

EPO

**Uniprot**

P07865

**Synonyms:**

Erythropoietin, EPO

**Immunogen:**

Recombinant Macaca fascicularis Erythropoietin protein (28-192AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A