

EPPK1 Antibody



PACO26521

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EPPK1

Uniprot

P58107

Synonyms:

Epiplakin (450 kDa epidermal antigen), EPPK1, EPIPL

Immunogen:

Recombinant Human Epiplakin protein (1-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A