

ESD Antibody, HRP conjugated



PACO26542

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serine hydrolase involved in the detoxification of formaldehyde.

Gene ID:

ESD

Uniprot

P10768

Synonyms:

S-formylglutathione hydrolase (FGH) (EC 3.1.2.12) (Esterase D) (Methylumbelliferyl-acetate deacetylase) (EC 3.1.1.56), ESD

Immunogen:

Recombinant Human S-formylglutathione hydrolase protein (2-282AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A