

Hemoglobin Antibody, FITC conjugated



PACO26593

Product Information

Size:

50ug

Reactivity:

Pigeon

Source:

Guinea pig

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

Synonyms:

Immunogen:

Pigeon Hemoglobin (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A