

FABP1 Antibody, FITC conjugated



PACO26623

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds free fatty acid, and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.

Gene ID:

FABP1

Uniprot

P07148

Synonyms:

Fatty acid, binding protein, liver (Fatty acid, binding protein 1) (Liver-type fatty acid, binding protein) (L-FABP), FABP1, FABPL

Immunogen:

Recombinant Human Fatty acid, binding protein, liver protein (1-127AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A