

# FABP1 Antibody, Biotin conjugated



PACO26624

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds free fatty acid, and their coenzyme A derivatives, bilirubin, and some other small molecules in the cytoplasm. May be involved in intracellular lipid transport.

**Gene ID:**

FABP1

**Uniprot**

P07148

**Synonyms:**

Fatty acid, binding protein, liver (Fatty acid, binding protein 1) (Liver-type fatty acid, binding protein) (L-FABP), FABP1, FABPL

**Immunogen:**

Recombinant Human Fatty acid, binding protein, liver protein (1-127AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A