

FABP3 Antibody, Biotin conjugated



PACO26628

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FABP3

Uniprot

P05413

Synonyms:

Fatty acid, binding protein, heart (Fatty acid, binding protein 3) (Heart-type fatty acid, binding protein) (H-FABP) (Mammary-derived growth inhibitor) (MDGI) (Muscle fatty acid, binding protein) (M-FABP), FABP3, FABP11 MDGI

Immunogen:

Recombinant Human Fatty acid, binding protein, heart protein (2-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A