

# ANXA6 Antibody, FITC conjugated



PACO26631

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May associate with CD21. May regulate the release of Ca<sup>2+</sup> from intracellular stores.

**Gene ID:**

ANXA6

**Uniprot**

P08133

**Synonyms:**

Annexin A6 (67 kDa calelectrin) (Annexin VI) (Annexin-6) (Calphobindin-II) (CPB-II) (Chromobindin-20) (Lipocortin VI) (Protein III) (p68) (p70), ANXA6, ANX6

**Immunogen:**

Recombinant Human Annexin A6 protein (2-245AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A