

ANXA6 Antibody, Biotin conjugated



PACO26632

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May associate with CD21. May regulate the release of Ca²⁺ from intracellular stores.

Gene ID:

ANXA6

Uniprot

P08133

Synonyms:

Annexin A6 (67 kDa calelectrin) (Annexin VI) (Annexin-6) (Calphobindin-II) (CPB-II) (Chromobindin-20) (Lipocortin VI) (Protein III) (p68) (p70), ANXA6, ANX6

Immunogen:

Recombinant Human Annexin A6 protein (2-245AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A