

NUTM2F Antibody



PACO26665

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NUTM2F

Uniprot

A1L443

Synonyms:

NUT family member 2F, NUTM2F, FAM22F

Immunogen:

Recombinant Human NUT family member 2F protein (450-690AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A