

# FBXO24 Antibody, HRP conjugated



PACO26702

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

**Gene ID:**

FBXO24

**Uniprot**

O75426

**Synonyms:**

F-box only protein 24, FBXO24, FBX24

**Immunogen:**

Recombinant Human F-box only protein 24 protein (1-200AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A