

FCSD2 Antibody, Biotin conjugated



PACO26716

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FCSD2

Uniprot

O94868

Synonyms:

F-BAR and double SH3 domains protein 2 (Carom) (SH3 multiple domains protein 3), FCSD2, KIAA0769 SH3MD3

Immunogen:

Recombinant Human F-BAR and double SH3 domains protein 2 protein (1-290AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A