

FFAR2 Antibody



PACO26725

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

FFAR2

Uniprot

O15552

Synonyms:

Free fatty acid, receptor 2 (G-protein coupled receptor 43), FFAR2, FFA2 GPCR43 GPR43


Immunogen:

Recombinant Human Free fatty acid, receptor 2 protein (277-330AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung tissue using PACO26725 at dilution of 1:100.



Immunofluorescent analysis of MCF-7 cells using PACO26725 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).