

FFAR2 Antibody, FITC conjugated



PACO26727

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	FFAR2
Source:	Uniprot
Rabbit	O15552
Isotype:	Synonyms:
IgG	Free fatty acid, receptor 2 (G-protein coupled receptor 43), FFAR2, FFA2 GPCR43 GPR43
Applications:	Immunogen:
ELISA	Recombinant Human Free fatty acid, receptor 2 protein (277-330AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A