

FGA Antibody



PACO26733

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

FGA

Uniprot

Q8WW76

Synonyms:

FGA protein, FGA


Immunogen:

Recombinant Human FGA protein (36-218AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry analysis of human kidney tissue using PACO26733 at dilution of 1:100.