

FGA Antibody, FITC conjugated



PACO26735

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FGA

Uniprot

Q8WW76

Synonyms:

FGA protein, FGA

Immunogen:

Recombinant Human FGA protein (36-218AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A