

# Fgf21 Antibody, FITC conjugated



PACO26751

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Stimulates glucose uptake in differentiated adipocytes via the induction of glucose transporter SLC2A1/GLUT1 expression (but not SLC2A4/GLUT4 expression). Activity probably requires the presence of KLB.

**Gene ID:**

Fgf21

**Uniprot**

Q9JJN1

**Synonyms:**

Fibroblast growth factor 21 (FGF-21), Fgf21

**Immunogen:**

Recombinant Mouse Fibroblast growth factor 21 protein (29-210AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A