

Fgf23 Antibody, HRP conjugated



PACO26754

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulator of phosphate homeostasis. Inhibits renal tubular phosphate transport by reducing SLC34A1 levels. Regulator of vitamin-D metabolism. Negatively regulates osteoblasts differentiation and. matrix mineralization. Acts directly on the parathyroid to decrease PTH secretion. Upregulates EGR1. expression in the presence of KL.

Gene ID:

Fgf23

Uniprot

Q8VI82

Synonyms:

Fibroblast growth factor 23 (FGF-23), Fgf23

Immunogen:

Recombinant Rat Fibroblast growth factor 23 protein (25-251AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A