

# Fgf23 Antibody, FITC conjugated



PACO26755

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Regulator of phosphate homeostasis. Inhibits renal tubular phosphate transport by reducing SLC34A1 levels. Regulator of vitamin-D metabolism. Negatively regulates osteoblasts differentiation and. matrix mineralization. Acts directly on the parathyroid to decrease PTH secretion. Upregulates EGR1. expression in the presence of KL.

**Gene ID:**

Fgf23

**Uniprot**

Q8VI82

**Synonyms:**

Fibroblast growth factor 23 (FGF-23), Fgf23

**Immunogen:**

Recombinant Rat Fibroblast growth factor 23 protein (25-251AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A